

soreFix®



Duo Patch

Cold sore patch with 0.5% hyaluronic acid

- Stimulates wound healing & relieves itching
- Reduces risk of infection
- Transparent, hypoallergenic, make-up resistant
- Effective for 12 hours

SoreFix Duo Patch is an ultra-thin, transparent, hypoallergenic patch which adheres to the lip. The patch contains hyaluronic acid which carefully moisturises the cold sore. This aids the recovery of the cold sore, relieves itching, creates a protective layer against external influences, reduces the risk of infection and prevents scarring. The patch is effective for about 12 hours, not visible and makeup can be applied over it.

What is herpes labialis?

Herpes labialis is a herpes simplex virus infection that appears as a blister on the lip. These filled blisters are called cold sores.

What are the symptoms of a herpes labialis infection?

The first symptoms are a warm, tingling sensation on the lip and redness. Fluid-filled blisters form within hours. A cold sore can be triggered by a weakened immune system, fever, menstruation, stress, exposure to intense sunlight, or a reaction to certain nutrients. A cold sore can always reoccur because the virus remains present and can be reactivated by the causes previously stated above.

What to avoid when you have herpes labialis?

If you have a cold sore, avoid lip contact with other people until the blisters have completely disappeared. Herpes labialis is highly contagious! Avoid very hot or very cold foods as these can make the symptoms worse.

How does SoreFix Duo Patch work?

The soft side of the patch is covered with a hypoallergenic adhesive layer in which the hyaluronic acid is evenly distributed. Hyaluronic acid has the property to store large amounts of water, creating a moist environment under the patch which stimulates wound healing for 12 hours and relieves itching. The patch also forms a barrier on the skin against external influences, reducing the risk of infection and protecting the skin from scarring.

**Usage (see illustrations):**

Wash hands before using the Patch and make sure the skin is clean and dry (the use of creams or oils hinders the adhesion of the patch to the skin).



Remove the patch from the protective paper (1, 2), remove the first half of the patch without touching the adhesive part of the patch with your hands (3, 4, 5).



Apply the patch to the affected area and remove the second half of the patch until it is fully adhered (6).



It is recommended to use the patch until the cold sore is completely healed.

**Warning:**

- Can cause an allergic reaction
- Keep out of reach of children

Contents: 15 active patches

In case of serious incidents, report to the manufacturer and competent authority
SoreFix offers products that treat cold sores and mouth ulcers to help prevent cold sore outbreaks.
Goodbye sore, hello smile.



For more information, please visit www.sorefix.com

**A product from Serrix B.V.:**

Herengracht 458
1017 CA Amsterdam
The Netherlands

**Distributor:**

(Name)
(Address)
(Telephone)
(Mail)

**Eurosirel SpA**

Viale Europa, 30
Cusago 20047
Milano – Italy
eurosirel@eurosirel.com



L214EG102

REF CFUA701420